FR

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EN

Medical device of class 1 according to MDR (EU) 2017/745.
These devices are intended to be handled by anyone able to read and understand the instructions for use. **Patient profile** : Any type of patient whose pathology requires the use of a pressure sore prevention and/or treatment aid. **Operator's profile for installation, implementation, adjustments and maintenance** : Patient, healthcare professionals, maintenance technicians, third party. There are no particular restrictions on the handling and use of the device. **INTENDED USE OF PRODUCT**
Positioning aids in bed and/or aids in the prevention and treatment of pressure sore/bedsores. **INDICATION**
Refer to each item indication. The choice of device must be approved by the healthcare professional. **CONDITIONS OF USE**
Devices with an integral POLYMAILLE® / POLYMAILLEHD can be used for single or multiple patients. They can be reconditioned in accordance with the product intructions and the cleaning and disinfection instructions. The first preventive measure against bedsores: the body should be raised as frequently as possible to limit compression of the blood vessels, thereby enabling a correct blood circulation.

CONTRA-INDICATIONS
Do not use the device beyond the intended use and indications. Do not use if you are allergic to polyurethane. For any problem, consult your healthcare provider.

SIDE EFFECTS
There are no undesirable side effects. However, the patient may experience the effects inherent to immersion: possible (limited) reduction in perception and mobility. Increased sensation of heat.

PERFORMANCE
Allows relief of support areas through immersion or discharge. Helps maintaining lying-down posture.

CLINICAL BENEFIT
Improves tissue oxygenation in the contact zone.

INFORMATION FOR HEALTHCARE PROFESSIONAL
Use positioning aids or pressure relief devices in patient with pressure sore/bedsore. Make sure the device is suitable for the patient's environment. Ensure continuity of treatment.

STORAGE, HANDLING, TRANSPORT
Preferably flat. Do not expose close to a significant heat source.

Imperatively in a dry and temperate place protected from heat and light in protective bag.

Handle the device with care to avoid damaging the cover.

Transport in a transport bag or in its box.

For any handling / Transport please follow these instructions.

PRECAUTIONS OF USE
Do not use without cover.

Do not place any foreign body between the patient and the equipment besides clothing.

Be careful not to damage the cover. Frequent degradations observed: cigarette burns, pet scratches, needle holes, cuts with scalpels or sharp objects, damage caused by the use of equipments with sharp edges.

We recommend obtaining an additional protective cover.

Do not use spare covers, accessories or components other than those supplied by SYSTAM®.

Refer to "Cleaning & Disinfection" instructions.

If you observe deterioration of the product (foam deformation or leakage) or bottoming effect, contact your provider.

Keep away from flammable sources.

Report to your caregiver as soon as possible any abnormal event such as pain or redness or whitening of support points.

Any serious incident that has occurred in connection with the device should be notified to the manufacturer and to the local competent authority.

TECHNICAL SPECIFICATIONS

Cushion materials:

P910L: Moulded viscoelastic polyurethane (PU) foam + high resilience moulded polyurethane (PU) foam.

P911L: High resilience moulded polyurethane foam (PU).

P912L, P913L, P914L, and P915L, P918L, P919L: moulded viscoelastic polyurethane foam (PU).

Covers material:

POLYMAILLE®: polyurethane coating on jersey fabric®.

The lifespan of these products is 5 years when used in strict accordance with the care and maintenance instructions. The service life can be longer if the product is maintained with care and considerably shortened by extreme or improper use. This does not constitute an additional guarantee on these products.

GUARANTEE

Two years against any manufacturing defaults, from the date of purchase by the customer. Small imperfections (small tears, marbling, surface defects, etc.) due to the manufacturing process in no way affect the quality and effectiveness of the product and cannot be a reason for replacing the product under warranty. The warranty does not cover misuse of the product, normal wears of the device or of its cover. In the event of default, please contact your provider, and provide them the product lot number.

SAFETY & ENVIRONMENT

Do not burn or dispose of in the environment. Respect the recycling channels in place in your country. Do not ingest. 0 % C.F.C. foam. Complies with fire standards EN 597-1 & 597-2.

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