(FR)

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(EN)

Medical device of class 1 according to MDR (EU) 2017/745.
Regulation (EU) 2017/745.
These devices are intended to be handled by anyone able to read and understand the instructions for use.

Patient profile: Any patient whose pathology requires the use of a pressure sore prevention and/or treatment aid. **Operator's profile for installation, implementation, adjustments and maintenance**: Patient, healthcare professionals, maintenance technicians, third party. There are no particular restrictions on the handling and use of the device.

INTENDED USE OF PRODUCT
Postural aid in sitting position. Equipment for use in chair with armrests.

INDICATION
P951D: Back support and Therapeutic Combi-Packs can be used to improve positioning in wheel chair, improve seating orientation, reduce forward sliding and relieve back pain. **P951F**: Designed for patients presenting partial or complete postural deficits in sitting position following nerve or muscle damage. The P951F is positioning device is used to position the upper limbs in a sitting position. The choice of device must be approved by the healthcare professional.

CONDITIONS OF USE
Devices with an integral POLYMAILLE® / POLYMAILLEHD cover can be used for single or multiple patients. They can be reconditioned in accordance with the product instructions and the cleaning and disinfection instructions. **CONTRA-INDICATIONS**
Do not use the device beyond the intended use and indications. Do not use if you are allergic to polyurethane. For any problem, consult your healthcare provider. The device is not appropriated for using in bed. **SIDE EFFECTS**
There are no undesirable side effects. However, the patient may experience the effects inherent to immersion: possible (limited) reduction in perception and mobility. Increased sensation of heat.

DEVEICE PERFORMANCES
Maintaining posture in sitting position. **CLINICAL BENEFIT**
P951D: Improved patient comfort. Contributes to the reduction of back pain and to the stability of the trunk. **P961F**: Reduces the risk of shoulder subluxation and contributes to the maintenance/stability of trunk balance.

INFORMATION FOR HEALTHCARE PROFESSIONAL
Make sure the device is suitable for the patient's environment. Ensure continuity of treatment. **STORAGE, HANDLING, TRANSPORT**
- Preferably flat. Do not expose close to a significant heat source. - Store in a protective bag, in a dry temperate place away from heat and light.

- Handle the device with care to avoid damaging the cover. - Transport in a transport bag or in its box.

PRECAUTIONS OF USE
- Do not use without cover. - Do not place any foreign body between the patient and the equipment besides clothing.

- Be careful not to damage the cover. Frequent degradations observed: scapulae burns, pet scratches, needle holes, cuts with scalpels or sharp objects, damage caused by the use of equipment with sharp edges. - We recommend obtaining an additional protective cover as a spare cover.

- Do not use spare covers, accessories or components other than those supplied by SYSTAM®. - Refer to «Cleaning & Disinfection» instructions.

- If you observe deterioration of the product (foam deformation or leakage) or bottoming effect, contact your provider.

- Keep away from flammable sources. - Report to your caregiver as soon as possible any abnormal event such as pain or redness or whitening of support points.

- Any serious incident that has occurred in connection with the device should be notified to the manufacturer and to the local competent authority.

TECHNICAL SPECIFICATIONS

P951D / P961F: Polyurethane (PU) viscoelastic foam with memory effect (80 kg/m³).

P952D: High resilience polyurethane (PU) foam (40 kg/m³).

Covers material:
POLYMAILLE®: Polyurethane coating on jersey fabric.

EXPECTED LIFETIME
The lifespan of these products is 5 years when used in strict accordance with the care and maintenance instructions. Product lifetime can be extended by careful maintenance and considerably shortened by extreme or inappropriate use. This does not constitute an additional guarantee on these products.

GUARANTEE
Two years against any manufacturing defaults, from the date of purchase by the customer. Small imperfections (small tears, marbling, surface defects, etc.) due to the manufacturing process in no way affect the quality and effectiveness of the product and cannot be a reason for replacing the product under warranty. The warranty does not cover misuse of the support, normal wear of the device or of its cover, in the event of default, please contact your provider, and provide the product lot number.

SAFETY & ENVIRONMENT
Do not burn or dispose of in the environment. Respect the recycling labels in place in your country. Do not ingest. 0 % C.F.C. foam. Complies with fire standards EN 597-1 & 597-2.

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